



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

POP Bandage Specifications

Catalogue No. /Reference No. /Products Code	JC04011
Material	Interlocked 100% cotton leno gauze, 40's, 20 threads, Gypsum Powder (calcium sulfate) and a plastic spool
Weight per Sqm	400~420GSM;
Color	White or other colors with non-toxic food pigment
Size	Width: 5cm, 6cm, 7.5cm, 8cm, 10cm, 12/12.5cm, 14cm, 15cm, 20cm, etc Length: 2m, 2.7m, 3m, 5yds, etc
Setting Time	2-4 minutes for Extra Fast Setting Bandages; 4-8 minutes for Fast Setting Bandages
Package	1 Roll/OPP bag
Sterile or Not	No Sterile
Rules & Classification	Rule 1, Class I





Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

POP Bandage Specifications

Reference No/ Product Code	Description	Size	Unit Package	Outside/Carton Package
JC040111	Plaster of Paris Bandage, 420gms, fast setting time, white color, 5/6/7.5/8/10/12/12.5/15/20cm x 2m/2.7m/3m/4.6m;	5CM X 2/2.7/3/4.6m	1 Roll/Opp Bag, 12 Roll/zipper Bag	240 Rolls
JC040112		6CM X 2/2.7/3/4.6m		240 Rolls
JC040113		7.5CM X 2/2.7/3/4.6m		240 Rolls
JC040114		8CM X 2/2.7/3/4.6m		240 Rolls
JC040115		10CM X 2/2.7/3/4.6m		120 Rolls
JC040116		12/12.5CM X 2/2.7/3/4.6m		120 Rolls
JC040117		15CM X 2/2.7/3/4.6m		120 Rolls
JC040118		20CM X 2/2.7/3/4.6m	1 Roll/Opp Bag, 6 Roll/zipper Bag	60 Rolls



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

POP Bandage Specifications

Reference No/ Product Code	Description	Size	Unit Package	Outside/Carton Package
JC040119	Plaster of Paris Bandage, 420gms, fast setting time, colored, 5/7.5/10/15/20cm x 2m/2.7m/3m/4.6m;	5CM X 2/2.7/3/4.6m	1 Roll/Opp Bag, 12 Roll/zipper Bag	240 Rolls
JC0401110		7.5CM X 2/2.7/3/4.6m		240 Rolls
JC0401111		10CM X 2/2.7/3/4.6m		120 Rolls
JC0401112		15CM X 2/2.7/3/4.6m		120 Rolls
JC0401113		20CM X 2/2.7/3/4.6m	1 Roll/Opp Bag, 6 Roll/zipper Bag	60 Rolls