



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,
Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Zinc Oxide Adhesive Plaster Specifications

Catalogue No. /Reference No. /Products Code	JC05001
Material	100% cotton fabric or cotton and polyester fabric coated with zinc oxide adhesive or hot-melt glue
Color	White, Skin Color, Yellow, Tan, etc
Size	Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm Length: 5m, 10m, 4yds, 5yds, 10yds
Package	1 Roll/Plastic (or Metal) Spool and Cover or Simple Package; 24/12/6 Rolls/Box; 30 Boxes/Carton
Shelf Life	3 Years
Sterile or Not	No Sterile
Rules & Classification	Rule 1, Class I





Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Zinc Oxide Adhesive Plaster Specifications

Reference No/ Product Code	Description	Size	Unit Package	Outside/Carton Package
JC050011	Zinc Oxide Adhesive Plaster, 1.25/2.5/5/7.5/10cm x 5m, with plastic spool and cover	White Color	24 Rolls/Box	720 Rolls
JC050012			12 Rolls/Box	360 Rolls
JC050013			6 Rolls/Box	180 Rolls
JC050014			6 Rolls/Box	180 Rolls
JC050015			6 Rolls/Box	180 Rolls
JC050016	Zinc Oxide Adhesive Plaster, 1.25/2.5/5/7.5/10cm x 5m/10yds, simple package, without plastic spool and cover	White Color	24 Rolls/Box	720 Rolls
JC050017			12 Rolls/Box	360 Rolls
JC050018			6 Rolls/Box	180 Rolls
JC050019			6 Rolls/Box	180 Rolls
JC0500110			6 Rolls/Box	180 Rolls



Shanghai Joy Crown Industry Co., Ltd.

Add: Rm.1105 No.55-56 Lane 21 Pu San Road, Pu Dong New Area,

Shanghai 200125, P.R. China

Tel: +86-21-58751019 Fax: +86-21-50907607 Mobile/Whatsapp: 008613818337572

Email: sales@joycrown.com <https://www.medicalsupplyment.com>

Zinc Oxide Adhesive Plaster Specifications

Reference No/ Product Code	Description	Size	Unit Package	Outside/Carton Package
JC0500111	Zinc Oxide Adhesive Plaster, 1.25/2.5/5/7/5/10cm x 5m, with plastic spool and cover	Skin Color	24 Rolls/Box	720 Rolls
JC0500112			12 Rolls/Box	360 Rolls
JC0500113			6 Rolls/Box	180 Rolls
JC0500114			6 Rolls/Box	180 Rolls
JC0500115			6 Rolls/Box	180 Rolls
JC0500116	Zinc Oxide Adhesive Plaster, 1.25/2.5/5/7/5/10cm x 5m/10yds, simple package, without plastic spool and cover	Skin Color	24 Rolls/Box	720 Rolls
JC0500117			12 Rolls/Box	360 Rolls
JC0500118			6 Rolls/Box	180 Rolls
JC0500119			6 Rolls/Box	180 Rolls
JC0500120			6 Rolls/Box	180 Rolls